



Toxicology Services
1111 Newton Street
Gretna, Louisiana 70053 USA
800-433-3823 Toll-free Phone
504-361-8989 Phone
504-361-8298 Fax

Configurations for ToxScreen Integrated Split Key Drug Test Cup

ToxScreen Split Key Cup – 3-drug (KC330)		
Abbreviation	Drug / Drug Metabolite	Cut-off Level
COC	Cocaine metabolite	300 ng/mL
MAMP	Methamphetamine	1000 ng/mL
THC	Marijuana metabolite	50 ng/mL

ToxScreen Split Key Cup – 5-drug (KC500)		
Abbreviation	Drug / Drug Metabolite	Cut-off Level
AMP	Amphetamines	1000 ng/mL
COC	Cocaine metabolite	300 ng/mL
OPI	Opiates	2000 ng/mL
PCP	Phencyclidine	25 ng/mL
THC	Marijuana metabolite	50 ng/mL

ToxScreen Split Key Cup – 5-drug (KC530)		
Abbreviation	Drug / Drug Metabolite	Cut-off Level
COC	Cocaine metabolite	300 ng/mL
MAMP	Methamphetamine	1000 ng/mL
OPI	Opiates	2000 ng/mL
PCP	Phencyclidine	25 ng/mL
THC	Marijuana metabolite	50 ng/mL

ToxScreen Split Key Cup – 5-drug (KC540)		
Abbreviation	Drug / Drug Metabolite	Cut-off Level
AMP	Amphetamines	1000 ng/mL
COC	Cocaine metabolite	300 ng/mL
MAMP	Methamphetamine	1000 ng/mL
OPI	Opiates	2000 ng/mL
THC	Marijuana metabolite	50 ng/mL

ToxScreen Split Key Cup – 5-drug with specimen validity (KC500A)		
Abbreviation	Drug / Drug Metabolite	Cut-off Level
AMP	Amphetamines	1000 ng/mL
COC	Cocaine metabolite	300 ng/mL
OPI	Opiates	2000 ng/mL
PCP	Phencyclidine	25 ng/mL
THC	Marijuana metabolite	50 ng/mL
Specimen Validity Testing (pH, SG, & OX)		

Configurations for ToxScreen Integrated Split Key Drug Test Cup

ToxScreen Split Key Cup – 5-drug with specimen validity (KC530A)		
Abbreviation	Drug / Drug Metabolite	Cut-off Level
COC	Cocaine metabolite	300 ng/mL
MAMP	Methamphetamine	1000 ng/mL
OPI	Opiates	2000 ng/mL
PCP	Phencyclidine	25 ng/mL
THC	Marijuana metabolite	50 ng/mL
Specimen Validity Testing (pH, SG, & OX)		

ToxScreen Split Key Cup – 5-drug with specimen validity (KC540A)		
Abbreviation	Drug / Drug Metabolite	Cut-off Level
AMP	Amphetamines	1000 ng/mL
COC	Cocaine metabolite	300 ng/mL
MAMP	Methamphetamine	1000 ng/mL
OPI	Opiates	2000 ng/mL
THC	Marijuana metabolite	50 ng/mL
Specimen Validity Testing (pH, SG, & OX)		

ToxScreen Split Key Cup – 7-drug (KC730)		
Abbreviation	Drug / Drug Metabolite	Cut-off Level
BAR	Barbiturates	300 ng/mL
BZO	Benzodiazepines	300 ng/mL
COC	Cocaine metabolite	300 ng/mL
MAMP	Methamphetamine	1000 ng/mL
OPI	Opiates	2000 ng/mL
PCP	Phencyclidine	25 ng/mL
THC	Marijuana metabolite	50 ng/mL

ToxScreen Split Key Cup – 7-drug with specimen validity (KC730A)		
Abbreviation	Drug / Drug Metabolite	Cut-off Level
BAR	Barbiturates	300 ng/mL
BZO	Benzodiazepines	300 ng/mL
COC	Cocaine metabolite	300 ng/mL
MAMP	Methamphetamine	1000 ng/mL
OPI	Opiates	2000 ng/mL
PCP	Phencyclidine	25 ng/mL
THC	Marijuana metabolite	50 ng/mL
Specimen Validity Testing (pH, SG, & OX)		

ToxScreen Split Key Cup – 9-drug (KC990)		
Abbreviation	Drug / Drug Metabolite	Cut-off Level
AMP	Amphetamines	1000 ng/mL
BZO	Benzodiazepines	300 ng/mL
COC	Cocaine metabolite	300 ng/mL
MAMP	Methamphetamine	1000 ng/mL
MDMA	Ecstasy	500 ng/mL
MOR	Opiates	300 ng/mL
MTD	Methadone	300 ng/mL
PCP	Phencyclidine	25 ng/mL
THC	Marijuana metabolite	50 ng/mL

Configurations for ToxScreen Integrated Split Key Drug Test Cup

ToxScreen Split Key Cup – 10-drug (KC1000)		
Abbreviation	Drug / Drug Metabolite	Cut-off Level
AMP	Amphetamines	1000 ng/mL
BAR	Barbiturates	300 ng/mL
BZO	Benzodiazepines	300 ng/mL
COC	Cocaine metabolite	300 ng/mL
MAMP	Methamphetamine	1000 ng/mL
MDMA	Ecstasy	500 ng/mL
OPI	Opiates	2000 ng/mL
MTD	Methadone	300 ng/mL
PCP	Phencyclidine	25 ng/mL
THC	Marijuana metabolite	50 ng/mL