

AlereAutoLink

Experience and innovation have given Alere Toxicology Services, Inc. a unique advantage in providing accurate, dependable results that meet our clients' time requirements.

Alere's reporting and information website, www.AlereDataLink.com, links you with your drug screen data and enables you to more efficiently and effectively manage your own substance abuse testing program. This website also has an automated component that enables users to retrieve their drug screen data without ever manually logging into the website.

Through a full-time connection to the Internet, this program, called AlereAutoLink, enables users to perform print and export functions automatically. Once installed on the user's computer, the program prints or downloads test results at predetermined intervals throughout the day.

AlereAutoLink is easily configurable so that you can get your results how and when you want them, and after the initial set-up, no additional action is required by the user. It is also an efficient and reliable substitute for computer fax. You can set the program to print our download results at intervals of your choosing, from every 10 minutes to once every 24 hours.

AlereAutoPrint

Users of AlereAutoPrint can print drug test results in one of two formats – list or report. The List format provides a listing of the drug tests that have been completed and includes limited information on each result. The Report format provides a detailed and easy to follow report on each test result that includes all the necessary information regarding the test, making it a certified

copy of the result. Once the user has selected a format and interval for printing, AlereAutoPrint will search Alere's website for completed results and print them in the format selected. When AlereAutoPrint has printed a result, it will mark the result, and at the next printing interval, that result will not be printed, eliminating duplicate reports.

AlereAutoExport

Those clients that require our results to be imported into a database or software system can utilize AlereAutoExport to set intervals throughout the day to export results from Alere. These results can be downloaded in a ASCII or XML format to the user's computer until he/she is ready to import them. After a user sets an interval for exporting, AlereAutoExport will search Alere's website for completed results, compile them into a single file, and save them to the client's computer. As with printing, duplicate records are eliminated by marking previously exported results and not exporting them again.

Requirements for AlereAutoLink

- Windows 98, NT 4.0, 2000 or XP
- High-speed Internet Connection
- Internet Explorer 5.01 or higher
- 50Mb of free disk space
- Pentium 166 or Higher