

Tox Screen Drugs of Abuse Test Card

Alere offers Point of Collection Testing (POCT) devices to detect drugs of abuse in human urine. Alere has offered POCT devices since the early 1990's and reviews its products each year to provide clients with items that live up to Alere's high standards.

The FDA-registered Tox Screen card is a one-step immunoassay test for the qualitative detection of drugs in human urine. Simply collect a urine specimen in a sterile cup, dip the card into the specimen, and read the results when they have developed.

While the test is negative specific, it only provides a preliminary, non-quantitative positive. A more specific, alternative method, such as Gas Chromatography/Mass Spectrometry (GC/MS), should be used to confirm on-site positive results. Should results of the card indicate a positive, Alere recommends that the collected specimen be forwarded to our laboratory for analysis.

Available Configurations:

Alere offers several configurations of the Tox Screen Card, including single-drug devices for most substances and combinations ranging from 2 to 12 drugs on the same card.

Features:

- Easy-to-use, one-step device requires no incubation period, specimen transfers, or washing steps.
- Rapid, accurate, and sensitive "dip and read" test procedure ensures reliable test results.
- Test controls are built into the device.
- Results appear clear and are easy-to-read.
- Results are available in as few as 3 minutes.
- 90 mL sterile cups with temperature strip are available to use with the device.
- Devices may be stored at room temperature for 18 months.
- Devices must be ordered in increments of 25.

