

Tox Screen Drugs of Abuse Testing Cup

Alere offers Point of Collection Testing (POCT) devices to detect drugs of abuse in human urine. The Tox Screen Cup is an integrated collection cup and testing device that does not require any handling or manipulation to activate the test.

The FDA-registered Tox Screen cup allows for easy, single-step urine collection and testing while still keeping the device tamper-proof. Each cup is self-contained for easy transfer to our laboratory for confirmatory testing and includes a temperature and identification label to ensure specimen integrity and documentation.

In addition to providing quality POCT devices, Alere is committed to client support. When you are sold a product from Alere, you can also receive information on the use of that product and access to our Client Services Department to assist you in understanding the testing process and interpreting test results.

Available Configurations:

Alere offers several configurations of the Tox Screen Drugs of Abuse Testing Cup, including combinations of three to seven drugs and specimen validity testing.



Features:

- Integrated collection cup and test device allows for easy urine collection and prevents tampering.
- Each cup has a temperature strip affixed to it to enable users to verify the validity of the specimen as the urine is provided.
- Test results are kept confidential until the authorized individual is ready to view them.
- Collectors may record pertinent information regarding the specimen directly on the collection cup to aid in identification.
- Test controls are built into the device.
- Easy-to-read results are available in 5 minutes.
- Devices may be stored at room temperature for 18 months.
- Devices must be ordered in increments of 25.