

Tox Screen Split Key Drug Test Cup

Alere offers Point of Collection Testing (POCT) devices to detect drugs of abuse in human urine. The Tox Screen Split Key Cup is an integrated collection cup and testing device with a separated sample well to split the urine for the device's testing.

The FDA-cleared Tox Screen Split Key Cup allows for easy, single-step urine collection and testing while still keeping the device tamper-proof. The device comes with a uniquely designed key that activates the testing procedure, so that the test can be performed at the collector's convenience.

When the key is inserted into the side of the device cup, a portion of the urine is separated to perform the analysis, leaving the remaining urine unaffected by the testing procedure. The testing device draws urine from the bottom of the split sample well, eliminating the need to tip or shake the device.

Available Configurations:

Alere offers several configurations of the Tox Screen Integrated Split Key Cup, including combinations of four to 10 drugs and specimen validity testing.



Features:

- Integrated collection cup and test cassette allows for easy urine collection and prevents tampering.
- Test results are kept confidential until the authorized individual is ready to view them.
- Split specimen for the testing procedure ensures remaining specimen is unaffected for confirmatory testing.
- Test devices are self-contained and portable and can be sent directly to our laboratory for confirmation.
- Devices may be stored at room temperature for 18 months.
- Test controls are built into the device.
- Easy-to-read results are available in 3 to 8 minutes.
- Devices must be ordered in increments of 25.